

EXPERIMENT NO - 2

OBJECT:- To determine the chloride content $[Cl^-]$ in the supplied water sample by Mohr's Method.

REQUIREMENT:- Burette, pipette, conical flask, measuring flask funnel, standard silver nitrate ($N/50$), distilled water.

INDICATORS:- Potassium chromate (K_2CrO_4) indicator.

PRINCIPLE:- This example is an example of precipitation titration. Mohr's method is used to determine chloride content in given sample. In this method a slightly alkaline solution solution is titrated against standard silver nitrate solution using potassium chromate as an indicator, as the titration proceeds, the chloride ions present in the sample react with $AgNO_3$, forming insoluble white precipitate of silver chloride. The extra drop of $AgNO_3$ white precipitate reacts with Indicator forming a brick red precipitate of silver chromate. The change of colour from yellow to brick red marks the end point.

OBSERVATION:- (1) Concordant reading = 0.2 ml

(2) Concordant reading = 4.5 ml



OBSERVATION :- (1) Reading with blank sample

(1)

S.NO	Volume of distilled water in (ml)	Burette reading (ml) initial	Burette reading (ml) final	Vol. of AgNO_3 used in (ml)
1	10	0.0	0.2	0.2
2	10	0.2	0.4	0.2
3	10	0.4	0.6	0.2

$$\text{concordant reading} = V_A = \underline{0.2 \text{ ml}}$$



(2) Reading with sample solution.

S.NO	Vol. of Sample solution taken (ml)	Burette Reading initial	Burette Reading final	Vol. of AgNO_3 used in (ml)
1	10	0	4.5	4.5
2	10	4.5	9	4.5
3	10	9	13.5	4.5

$$\text{Concordant reading} = V_B = \underline{4.5 \text{ ml}}$$

SOLUBLE NITRATES

REACTION :-



CALCULATIONS:-

Volume of sample solution = 10 ml.

$$\text{Vol. of } \text{AgNO}_3 \text{ used} = V_B - V_A = V_{\text{ml}} = 4.5 - 0.2 = 4.3 \text{ ml}$$

$$\text{We know that, } N_1 V_1 = N_2 V_2$$

(AgNO_3) Sample solution.

$$\frac{1}{50} \times 4.3 = N_2 \times 10$$

$$\text{so, } N_2 = \frac{4.3}{500} = 0.0086$$

The strength of chloride ions in g/L =

Normality \times equivalent weight of Cl.

$$\Rightarrow \frac{4.3}{500} \times 35.5 \text{ gm/l} = 0.3053 \text{ g/l}$$

$$\Rightarrow \frac{4.3 \times 35.5 \times 10^3}{500} \text{ mg/l} = 305.3 \text{ mg/l}$$

RESULT :- The chloride content in given water sample is 305.3 ppm

PROCEDURE :- Take 10 ml of distilled water in a conical flask and add atleast 5 drop of K_2CrO_4 indicator. slowly add standard solution of AgNO_3 from the burette and the volume of AgNO_3 is noted at the end point. ie. when yellow colour starts changing to red coloured precipitate. Repeat the titration for concordant reading. Now take 10ml

of the given water sample in a conical flask. Add five drops of freshly prepared K_2CrO_4 solution. Titrate against $AgNO_3$ till permanganate brick red colour is obtained. Repeat the titration till concordant reading is obtained.

PRECAUTIONS:-

- (1) The apparatus should be washed with distilled water.
- (2) The reaction mixture should be properly shaken during titration.
- (3) Volume of Indicator should be same in all titration.